## TECHNISCHE INFORMATIONEN

DE

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EN

## 1. INTENDED USE

The 101-PL Pin Lock Unit is designed for a fixing of the pin of a silicone or gel socket with thread M10 within the prosthesis. The Pin Lock Unit consists of a mechanism with worm gear. This worm gear allows locking the locking pin in any position and to pull-in the socket into the prosthesis by turning on the outside release button. This mechanism of the Pin Lock Unit lies open on the outside of the prosthetic socket to enable an easy access at any time. To connect to the modular system, the pin-lock unit has a 4-hole Euro connection. The use of this 101-PL Pin Lock Unit is limited to a body weight of 150 kg. Additional loads and patient activity are also to be considered for weight selection.

## 2. PACKAGING CONTENT

The 101-PL Pin Lock Unit consists of:

- Clutch Lock Locking Device
- 4 x Countersunk Bolt DIN 7991 M6x40 12.9
- ca. 1,2 rmt 901-MT Molding Textile

## 3. ASSEMBLY

Only use the provided bolts, tighten them to a torque of 12Nm and secure the threads.

## 4. RISKS

Point out to the patient that correctly aligning of the locking pin inside the locking device is very important. The pin should not be inserted into the mechanism with full strength. If the pin is not correctly aligned to the mechanism, there is the risk that the pin could become stuck in the mechanism and may be very difficult to loosen. The functioning of the mechanism is to be checked in regular intervals, at least once a year.

## 4. WARRANTY

This warranty applies to faults in material and workmanship and is granted to orthopedic workshops only if all given proceeding instructions are complied with and the products are used if intended. The functioning of the mechanism is to be checked in regular intervals but at least once a year. The product has been tested according to ISO recommendations. Amparo GmbH offers a warranty of two (2) years in case of appropriate use.

## 5. DECLARATION OF CONFORMITY

We declare that the 101-PL Pin Lock Unit complies with the requirements of the Guideline MDR 2017/745.

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